



Acme Sanitary Ware Co Ltd

Acme Building, 22-28 Nanking Street, Kowloon, Hong Kong
 Tel: 2388 7171 Fax: 2710 8012
 E-mail: acme@acmesanitary.com.hk
 Website: www.acmesanitaryware.com.hk

PROJECT		REF		REV	ITEM CODE	
LOCATION		DATE			PAGE	

SANITARY WARE SPECIFICATION SHEET

Item Descriptions	Toilet paper dispenser, standard roll, steel white epoxy	<div>Illustration/ Drawing</div> <div><p>* Technical Specifications</p><table><tr><td>Toilet paper roll capacity</td><td>3 u</td></tr><tr><td>Dimensions</td><td>380x 116x 130 mm</td></tr><tr><td>Weight (empty dispenser)</td><td>1,3 Kg</td></tr><tr><td>Toilet paper roll diameter</td><td>Ø 105-125 mm (maximum)</td></tr><tr><td>Size toilet paper roll</td><td>110 mm (maximum)</td></tr><tr><td>Recommended installation height to the bottom</td><td>65-75 cm</td></tr></table></div>	Toilet paper roll capacity	3 u	Dimensions	380x 116x 130 mm	Weight (empty dispenser)	1,3 Kg	Toilet paper roll diameter	Ø 105-125 mm (maximum)	Size toilet paper roll	110 mm (maximum)	Recommended installation height to the bottom	65-75 cm
Toilet paper roll capacity	3 u													
Dimensions	380x 116x 130 mm													
Weight (empty dispenser)	1,3 Kg													
Toilet paper roll diameter	Ø 105-125 mm (maximum)													
Size toilet paper roll	110 mm (maximum)													
Recommended installation height to the bottom	65-75 cm													
Dimensions	L116 x W130 x H380 mm													
Model	PR0781													
Material/ Finish	Stainless Steel/ White epoxy													
Manufacturer	Mediclinics, S.A. (Spain)													
Supplier	Acme Sanitary Ware Co. Ltd Mr. Eric Wong/ Mr. Wilson Hung													
Contact Tel/Fax	(852) 2388-7171 / (852) 2710-8012													
E-mail	acme@acmesanitary.com.hk													
Website	www.acmesanitary.com.hk													
<div>Components and Materials</div> <ul style="list-style-type: none">Steel white epoxy coated one-piece body, 0,8 mm thick, fully sealed. The supports of the screws in the back of the appliance, are acting as an anti-moisture support, as well as slots to avoid the screws damaging the toilet paper. It is mounted on the wall by means of 3 wood screws with their corresponding 5 mm wall plugs. It includes a lock system with standard mediclinics ® key, that allows opening the lid for replenishment.Steel white epoxy coated one-piece seamless lid, 0,8 mm thick, fully sealed. Fixed to the body by means of 2 rivets that allows swinging of the lid for the replenishment through the top opening. The extra tissue rolls descend when the previous one has been finished, by a simple hand manipulation through the side window. There is no need to open the lid for the removal of consumed rolls.Frontal sight-hole and slots on each side of body indicate refill time. <div>Operation</div> <p>Pull visible toilet paper, at your own convenience.</p>														

Note: